


2203R

Flexible Film Tape



Backing: 3 mil clear polyethylene film with matte

finish

Adhesive: Non-sensitizing medical-grade acrylic

(M109)

Liner: 60lb kraft, silicone release treated on

one side

Applications

Device overlay

Medical device attachment

• Incise drape for surgical incision

Drape to patient attachment

Features Benefits

PE Film Backing	Matte finish for minimal glare	
	Highly conformable	
Medical-grade Acrylic Adhesive	 Excellent adhesion to skin and drape materials 	
60lb Release Liner	 Easy to die-cut and convert 	

Typical Test Values

	Standard UOM	Metric UOM	Test Method		
Total Thickness w/liner	7.7 mils	0.196 mm	PSTC 133		
Adhesive Thickness	1.2 mils	0.030 mm	PSTC 133		
Carrier Thickness	3.0 mils	0.076 mm	PSTC 133		
Liner Thickness	3.5 mils	0.089 mm	PSTC 133		
Adhesion to Steel	140 oz/in	15.32 N/cm	PSTC 101 Mod		
Product Safety					

The base adhesives in this medical product have been tested in accordance with ISO 10993 and found to be non-cytotoxic, non-irritating, and non-sensitizing.

In accordance with FDA cGMP requirements, Berry reminds the customer that ultimate fitness for use must be determined by the manufacturer of the finished device.

Recommended Storage Conditions: 40-60% Humidity, 60 °F - 80 °F

The information presented herein was prepared at BERRY GLOBAL, INC. (25 Forge Parkway Franklin, MA 02038) by qualified technical personnel. To BERRY's knowledge it is true and accurate. However, the information and recommendations are furnished for these products with the understanding that the purchaser or user, as the case may be, will independently determine that the product is suitable for the intended use, and that such use complies with all applicable federal, state and local laws and regulations. The data are submitted only for the user's information and consideration, and do not constitute a warranty of any kind (including but not limited to a warranty of merchantability or of fitness for a particular purpose), or a representation for which BERRY



